

## PRODUCT INFORMATION

**Product Name :** DynaMarker RNA Low II Easy Load  
**Code No. :** DM157  
**Range :** 20-500 base of RNA  
**Size :** 125 µl, about 25 loadings  
**Loading :** 5 µl is recommended for loading to a well  
(0.1 µg of each RNA / 5 µl )

### Description :

The DynaMarker RNA Low II Easy Load is supplied in a ready-to-use mixture of loading dye (containing formamide, EDTA sodium salt, bromphenol blue) and RNAs. It is prepared for denaturing polyacrylamide gel electrophoresis but not agarose gel electrophoresis. The DynaMarker RNA Low II Easy Load has seven single-stranded RNAs, 20, 50, 100, 200, 300, 400 and 500 bases. The 20-base and 50-base RNA are synthesized by chemically (not phosphorylated), others are synthesized by *in vitro* transcription. In 5 µl of the DynaMarker RNA Low II, each RNA amount is approximately 100 ng. It is useful for estimating RNA amount approximately. The DynaMarker RNA Low II Easy Load can be visualized by UV light after ethidium bromide staining.

### Storage condition :

Store at -80 °C.  
Repeated freeze/thaw cycles should be avoided.

### Quality Control :

After 18 hr incubation of the DynaMarker RNA Low II Easy Load at 37 °C, no visible degradation of the marker is observed in 5 % polyacrylamide / 8M urea gel electrophoresis.

### Supplied product : RNA Loading buffer PA


RNA Loading buffer PA is manufactured for denaturing polyacrylamide gel electrophoresis but not agarose gel electrophoresis. The loading buffer has a composition of 80 % formamide, 10 mM EDTA sodium salt (pH8.0), 0.025 % bromphenol blue. Store RNA Loading buffer PA at -80 °C. Repeated freeze/thaw cycles should be avoided. It is 1 × to 2 × solution. Use more than one volume of RNA solution.

### Note :

RNA is very sensitive to degradation by nucleases. To avoid damaging the DynaMarker RNA Low II Easy Load, use extreme care during manipulations to prevent nuclease contamination. Wear gloves and use clean apparatus. Glassware should be pretreated with diethyl pyrocarbonate (DEPC). Nuclease-free disposable plasticware should be used. Solutions and reagents to mix the product should be high grade and nuclease-free. To use, thaw the DynaMarker RNA Low II Easy Load on ice and keep it on ice while using. For heat denaturation, transfer aliquot of the DynaMarker RNA Low II Easy Load to another tube, then heat it. Avoid repeated heat denaturizing.

‡ Formamide is suspected to be harmful. It is irritate to the eyes and skin. Wear appropriate gloves and safety glasses. Put a lid tightly at the time of storage.

*This product is research use only*



**DynaMarker RNA Low II Easy Load**

Electrophoresis profile of  
DynaMarker RNA Low II Easy Load  
(5 µl) on 5 % of acrylamide, 8 M urea  
gel with 1 × TBE buffer as running  
buffer

## PRODUCT INFORMATION

---

### Recommended usage :

The <sup>DynaMarker</sup> RNA Low II Easy Load is manufactured for denaturing polyacrylamide gel electrophoresis. As recommended usage, <sup>DynaMarker</sup> RNA Low II Easy Load is run on 5 % polyacrylamide / 8M urea gel as below. Effective range of separation of RNAs is about 50 – 500 base in 5 % polyacrylamide / 8M urea gel.

### Procedure

#### 1. Preparation of 40 % Acrylamide : bis solution

Acrylamide	190 g
N, N-methylenebisacrylamide	10 g
H <sub>2</sub> O	to 500 ml

After mixing, filter the solution through a nitrocellulose filter (0.45 µm pore size).

#### 2. Preparation of 5 % polyacrylamide / 8M urea gel (20 ml gel)

40 % acrylamide : bis solution	2.5 ml
Urea	9.6 g
10 × TBE	2.0 ml
H <sub>2</sub> O	to 20 ml

After urea is dissolved completely, add 20 µl of TEMED and 160 µl of 10 % ammonium persulfate. Mix quickly and then pour the gel into the mold of a vertical gel apparatus (7 cm × 8 cm, thickness 1.0 mm). The gel apparatus should be assembled according to the manufacture's protocol and ready to run with 1 × TBE buffer.

#### 3. Loading and electrophoresis

Mix RNA to be analyzed (for example, RNA transcript) and RNA Loading buffer PA as below.

RNA sample	dried precipitate or 2 µl (0.5 – 2 µg)
<u>RNA Loading buffer PA</u>	<u>5 µl - - - over one volume of RNA sample</u>

Mix in a small tube, total 5-7 µl

Transfer aliquot (5-10 µl) of <sup>DynaMarker</sup> RNA Low II Easy Load to a small tube. Heat RNA mixed with RNA Loading buffer PA and <sup>DynaMarker</sup> RNA Low II Easy Load at 80 °C for 3 min, and transfer the tube on ice immediately, then load onto a well of 5 % polyacrylamide / 8M urea gel and start electrophoresis. After the tracking dye has migrated an appropriate distance through gel, stop the electrophoresis. To stain with ethidium bromide, disassemble the apparatus and transfer the polyacrylamide gel to a gel tray filled with 1 × TBE buffer containing 10 µg/ml ethidium bromide. Stained RNA can be visualized using UV transilluminator.

### Reference:

Sambrook, J. and Russell, D.W. (2001) Molecular Cloning: A Laboratory Manual, 3rd ed., Cold Spring Harbor Laboratory Press, Cold Spring Harbor, NY.