

PRODUCT INFORMATION

Product Name : 6 X BPB loading Dye

Code No.	Size
DM210	1 ml X 4
DM212	1 ml X 2

This product is research use only.

Description : 6 X BPB loading Dye is useful for loading DNA samples into wells on agarose gel electrophoresis. It contains EDTA• 2Na, glycerol, bromphenol blue and 10 mM Tris-HCl (pH8.0).

Storage condition : Stable at room temperature.
For long time storage, -20°C is recommended.

Recommended usage: Mix 6 X BPB loading Dye with 1/5 volume of DNA solution to a well on agarose gel electrophoresis.

DNA sample	5 μl
<u>6 X BPB loading Dye</u>	<u>1 μl</u> - - - addition of 1/5 volume of DNA sample
Mix	



Load to a well.