

PRODUCT INFORMATION

Product Name :	Sonicated Salmon Sperm DNA		
Code No. :	F012, F013		
Size :	1 ml / tube	F012	1 ml
		F013	5 ml (1ml × 5)
Concentration:	10 mg/ml		
Description:	Salmon Sperm DNA was shared by sonication. DNA is dissolved in 10 mM Tris-HCl (pH 8.0), 1 mM EDTA.		
Quality Control:	Intrinsic DNase activity was not detected in the product, by the use of pUC18 plasmid DNA as a control.		
Preparation:	Boiling for 3 minutes and cooling in an ice bath before use.		
Storage Condition:	- 20 °C		

For Research Use Only. Not for use in diagnostic procedures.

Ver. 4.0