


IN3-VEC31 pFAST-G01



Name	Position
Left border (LB)	1-26
NOS terminator (NOST)	117-370
HPT	568-1590
NOS promoter (PNOS)	1605-1875
OLE1 promoter (POLE1)	1920-3363
OLE1-GFP fusion protein (OLE1-GFP)	3364-4893
35S terminator (T35S)	5177-5294
attR1	5371-5495
CmR	5604-6263
ccdB	6605-6910
attR2	6951-7075

Name	Position
Right border (RB)	7185-7209