

PURE AMP Compatibility Chart	MTCPCRP_ P9602-N	MTCPCRP_ P9602-SRR	MTCPCRP_ P9602-S	MTCLPP_ P9601-SRR	MTCLPP_ P9601-S	MTCLPP_ P9601-SW	MTCLPT_ P3801-QF	MTCLPT_ P3801-QW	MTCLPT_ P3801-Q	MTPCRT_ P3802-Q	MTPCRT_ P3802-QW
	Product page GO!			Product page GO!			Product page GO!			Product page GO!	
BENCHMARK SCIENTIFIC											
TC9639	v		v							v	v
TC-32 Min										v	v
APPLIED BIOSYSTEMS											
7000		v								v	v
7300		v								v	v
7500		v								v	v
7500FAST				v			v	v	v		
7900HT Fast Model (0.1 ml)				v							
7900HT Standard (0.2ml)		v									
StepOne							v	v	v		
StepOne Plus				v			v	v	v		
ViiA7 Fast Model (0.1 ml)				v			v	v	v		
ViiA7 (0.2ml)		v								v	v
QuantStudio Fast Model (0.1 ml)				v			v	v	v		
QuantStudio (0.2ml)		v								v	v
QuantStudio 3/5 Fast (0.1 ml)				v			v	v	v		
QuantStudio 3/5 (0.2ml)		v								v	v
5700/7700		v	v								
3100/3130 xl (0.1 ml Fast)				v							
310/3100/3130 xl (0.2ml Standard)		v	v								
3700/3730 xl (0.1 ml Fast)				v							
3700/3730 xl (0.2ml Standard)		v	v								
3500(xl) (0.1 ml Fast 96 Well)				v			v	v	v		
3500(xl) (0.2ml Standard 96 Well)		v	v							v	v
ROCHE											
LightCycler 480 and 48011						v					
LightCycler 96						v	v	v	v		
LightCycler Nano									v		
BIORAD											
iCycler iQ	v		v							v	v
iQ5	v		v							v	v
MyiQ	v		v							v	v
Chromo4	v		v							v	v
Opticon I & II					v	v	v	v	v		
mini Opticon											
CFX96					v	v	v	v	v		
OTHER MANUFACTURERS											
StrataGene MX4000/3000P/3005P	v		v							v	v
Agilent AriaMx					v	v	v	v	v		
IT-IS Life Science MyGo Pro									v		
OpenqPCR									v		
Eppendorf MasterCycle ep RealPlex	v		v							v	v
Thermo PikoReal 24 Real-Time PCR							v	v	v		
Biometra Optical ThermalCycler	v		v							v	v



DIAGNOCINE

Call us: 201-681-7263 // support@diagnocine.com